



## Disposable Surgical Gown SMS PE Laminated Reinforced, Non Sterile

 FDA CE  
MADE  
IN TURKEY

### Product Details

**Product Description :** SMS PE Laminated Reinforced Surgical Gown

**Brand Name:** Dexxon Medical

**Model Nr:** DXNMD- DNSLRSG-09

**Certificated :** CE and ISO13485

**Certificate Nr:** NVA-EC-21022307

**Product Name:** Disposable SMS Laminated Reinforced Sterile Surgical Gown

**Material quality:** PE (polyethylene) film-coated non-woven fabric

**Type :** Laminated Reinforced Surgical Gown Level 3 AAMI

**Usage:** Reinforced Surgical Gown

**Sterilization:** Non-sterile

**Weight:** 45 -55 g / m<sup>2</sup>

**Classification:** Class I

**Regular Color:** Medical Blue

### Basic Info

**Name:** SMS PE Laminated Reinforced Surgical Gown

**Color :** Blue

**Size :** S.M.L.XL.XXL.XXXL

**Regular Material:** PP/SMMS



## Disposable Surgical Gown SMS PE Laminated Reinforced, Non Sterile

  **MADE  
IN TURKEY**

**Performance indicators:** Technical requirements for disposable polyethylene laminated reinforced isolation clothing products / EN 13795-1:2019; AAMI PB70: 2012

**Level :** Level 3

### **Feature:**

- SMS Non-Woven Fabric
- Excellent water repellence to protect surgeons from contamination.
- Prevention of cross contamination by reducing lint
- Resistant to liquids with its barrier feature
- Resistant to moisture and bacterial penetration
- Puncture, tear and abrasion resistant.
- Free from toxic substances
- Non-Sterile

### **Application:**

This gown can be used during the surgery operations.  
This product is disposable and it prevents the bacterial transfer between the medical team and the patient, reduces the occurring of possible contamination and infection between the patient and the medical team.

### **Precautions:**

The product should be stored in a clean, dry place with proper ventilation and out of direct sunlight at 0°C-30°C (its relative temperature should not exceed 75%).